



**SRM VALLIAMMAI ENGINEERING COLLEGE**

(An Autonomous Institution)  
SRM Nagar, Kattankulathur – 603 203



**DEPARTMENT OF CHEMISTRY**

**GE3121-Chemistry Laboratory Manual**  
**Regulations 2023**  
**Semester I**  
**(2024-25)**

Department of Chemistry, SRMVEC

Published by the Principal, SRM Valliammai Engineering College, SRM Nagar, Kattankulathur-603203.

Copyright @2024, by The Principal, SRM Valliammai Engineering College.

No part of this laboratory manual may be reproduced or distributed in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise or stored in a database or retrieval system without the publisher's prior written permission. The program listings (if any) may be entered, stored and executed in a computer system, but they may not be reproduced for publication.

**GENERAL INSTRUCTIONAL POLICIES TO BE FOLLOWED IN  
THE CHEMISTRY LABORATORY**

1. Students must wear a knee length (41-42 inch) laboratory coat while entering into the laboratory.
2. Students must wear closed leather shoe (black or brown colour) for safety purpose.
3. Students must always come to the chemistry lab with lab manual, observation note, and record note.
4. Students should bring waste cloth for cleaning their hands.
5. Students may wear chemical safety goggles, if it is required.
6. Students should report to the concerned chemistry lab as per the time-table schedule.
7. Students, who turn up late to the lab, will not be permitted to do the experiment.
8. Students are required to prepare thoroughly to perform the experiment before coming to chemistry laboratory
9. After completion of the experiment, students must get the signature of the concerned staff in-charge in the observation note.
10. The record note must be submitted to get signature of the previous week experiment to the concerned staff while coming to lab.
11. When the experiment is completed, students should clean the apparatus carefully and should return all the reagent bottles/components/instruments taken for the experiment.
12. Any damages to apparatus occurred during the experiment, should be brought to the notice of lab in-charge. The cost of the repair or new apparatus should be brought by the students.
13. Students should be present in the lab till the completion of allotted hours of lab.
14. Students should not sit on the work tables or keep hands for rest on tables for their own safety.

15. All doors and windows must be kept open while working in the laboratory. Fans and lights should be switched off, if they are not required.
16. Long hair should be tied when in the laboratory so that it will not catch on fire or come into contact with chemicals.
17. Do not enter the laboratory until your lab instructor is present.
18. While leaving the lab ensure that all gas tap connections, water taps, electrical switches such as fans, lights shall be kept in off mode. All instruments must be properly shut down as per instructions of the teacher in charge.
19. Do not eat, drink and chew gum or smoke in the laboratory.
20. Do not perform unauthorized experiments. If you see someone else doing something you think may be dangerous, tell him or her to stop and/or report the incident to your lab instructor.
21. In case of emergency, do not get panicky and be brave enough to handle the situation. Medical attention should be requested immediately. If skin comes into contact with acid solutions, it should be washed with water and sodium bicarbonate solution alternately and continuously till skin is smoothened. Similarly exposure to caustic alkalis may be treated with dilute acetic or boric acid solutions and water. Fire scar can be washed thoroughly with water in a tap in the sink and a little oil or pipette be applied. If eyes come into contact with chemicals, wash them with plenty of water and clean them with a dry fresh cloth.
22. Fire buckets containing sand and water or fire extinguishers should always be available in the lab in easily accessible places.

### **Vision of the Chemistry Department**

Bridge the engineering professionals in recent advancements in chemistry for the developments of new devices and materials for the lives of individuals becoming different and complex global society.

### **Mission of the Chemistry Department**

To draw up outstanding educators, employers and researchers and to make progress in the profession of education through teaching and learning and art of research on chemical science towards engineering & technology.

**CHEMISTRY LABORATORY: (Any five experiments to be conducted)****OBJECTIVES**

- To inculcate experimental skills to test basic understanding of water quality parameters, such as acidity, alkalinity, hardness, DO, TDS, chloride, and chlorine.
- To familiarise the students with electroanalytical techniques like conductometry and flame photometry to determine the impurities in aqueous solution.
- To find the various characteristics of domestic water.
- To understand the Pseudo first-order kinetics reaction.
- To make the students know the synthesis of nanoparticles.

1. Determination of chloride content of water sample by Argentometric method.
2. Determination of total, temporary & permanent hardness of water by EDTA method.
3. Determination of DO content of water sample by Winkler's method.
4. Preparation of  $\text{Na}_2\text{CO}_3$  as a primary standard and estimation of acidity of a water sample using the primary standard.
5. Determination of types and amount of alkalinity in water samples.
6. Estimation of available chlorine in bleaching powder solution.
7. Conductometric titration of barium chloride against sodium sulfate (Precipitation titration).
8. Estimation of sodium /potassium present in water using a flame photometer.
9. Estimation of TDS of a water sample by gravimetry.
10. Preparation of nanoparticles (Ag/Au/ $\text{TiO}_2$ / $\text{ZnO}$ / $\text{CuO}$ ).
11. Pseudo first-order kinetics- ester hydrolysis.

**TOTAL: 30 PERIODS****OUTCOMES****At the end of the course, the student should be able:**

1. To infer the quality of water samples for alkalinity, hardness, DO, TDS, chloride, and chlorine.
2. To apply the knowledge on the estimation of metal ions, acidity and its precipitation nature towards their process.
3. To recognize the threshold limit for various characteristics of domestic water.
4. To identify the simple method of synthesis of nanoparticles.
5. To understand the pseudo first-order kinetics reaction from ester hydrolysis.

**TEXTBOOKS**

1. Vogel's Textbook of Quantitative Chemical Analysis (8<sup>th</sup> Edition, 2014).
2. Suchi Tiwari, Engineering Chemistry Lab Manual, Scitech Publications (India) Pvt. Ltd. (2<sup>nd</sup> Edition, 2013).
3. Pushpendra Kumar, Laboratory Manual for Engineering Chemistry, Reyansh Authortopic Pvt. Ltd., (1<sup>st</sup> Edition, 2022).

## CONTENTS

<b>EX. NO.</b>	<b>NAME OF THE EXPERIMENTS</b>	<b>PAGE NO.</b>
1	Determination of chloride content of water sample by argentometric method (Mohr's Method)	3
2	Determination of total, temporary & permanent hardness of water by EDTA method	7
3	Determination of types and amount of alkalinity in water samples.	12
4	Conductometric titration of barium chloride against sodium sulfate (Precipitation titration).	18
5	Preparation of Na <sub>2</sub> CO <sub>3</sub> as a primary standard and estimation of acidity of a water sample using the primary standard.	20
6	Estimation of TDS of a water sample by gravimetry.	25
<b>EXPERIMENT – BEYOND THE SYLLABUS</b>		
7	Preparation of nanoparticles (AgNPs).	27

## Titration 1: Standardization of silver nitrate

### Standard NaCl vs AgNO<sub>3</sub>

Sl.No	Volume of std. NaCl (ml)	Burette reading (ml)		Volume of AgNO <sub>3</sub> (ml)	Concordant value (ml)	Indicator
		Initial	Final			
1.	20	0				K <sub>2</sub> CrO <sub>4</sub>
2.	20	0				

### Calculation of the strength of silver nitrate

Volume of standard NaCl

$$(V_2) = 20 \text{ ml}$$

Strength of standard NaCl

$$(N_2) = 0.01 \text{ N}$$

Volume of unknown AgNO<sub>3</sub>

$$(V_1) = \text{_____ ml}$$

Strength of unknown AgNO<sub>3</sub>

$$(N_1) = ?$$

According to volumetric formula,

$$V_1 N_1 = V_2 N_2$$

$$N_1 = V_2 \times N_2 / V_1$$

$$N_1 = 20 \text{ ml} \times 0.01 \text{ N} / V_1$$

Strength of AgNO<sub>3</sub> (N<sub>1</sub>)

$$= \text{_____ N}$$



# 1. DETERMINATION OF CHLORIDE CONTENT OF WATER SAMPLE BY ARGENTOMETRIC METHOD (MOHR'S METHOD)

Expt. No.	
-----------	--

Date			
------	--	--	--

## AIM


To estimate the amount of chloride ion present in the water sample by Argentometric method (Mohr's method). You are provided with standard NaCl solution of strength 0.01 N and a link solution of AgNO<sub>3</sub> (approximately 0.01N).

## CHEMICALS REQUIRED


Standard NaCl solution, AgNO<sub>3</sub> solution, Potassium chromate indicator.

## PRINCIPLE

Natural water contains chloride ions in the form of NaCl, KCl, CaCl<sub>2</sub>, MgCl<sub>2</sub>. The concentration of chloride ions more than 250 ppm is not desirable for drinking purposes. This determination is based on precipitation titration. When AgNO<sub>3</sub> solution is added to the water sample, in presence of K<sub>2</sub>CrO<sub>4</sub>, the chlorides present in it are precipitated first as AgCl.



When all the Cl<sup>-</sup> ions are removed, AgNO<sub>3</sub> added from the burette will react with K<sub>2</sub>CrO<sub>4</sub> to give a reddish brown colour due to silver chromate (Ag<sub>2</sub>CrO<sub>4</sub>). This is the end point.



## PROCEDURE

### Titration 1: Standardisation of AgNO<sub>3</sub>

The burette is washed well with distilled water and rinsed with small amount of AgNO<sub>3</sub> solution. The pipette is washed with distilled water and rinsed with small amount of standard NaCl solution. 20ml of this solution is pipetted out into a clean conical flask. 1 ml or 3 drops of K<sub>2</sub>CrO<sub>4</sub> indicator solution is added and titrated against AgNO<sub>3</sub> solution taken in the burette. The end point is the change of colour from yellow to reddish brown. The titration is repeated for concordant values.

## Titration 2: Estimation of Chloride

### Water sample vs Standard AgNO<sub>3</sub>

Sl.No	Volume of std. water sample (ml)	Burette reading (ml)		Volume of AgNO <sub>3</sub> (ml)	Concordant value (ml)	Indicator
		Initial	Final			
1.	20	0				K <sub>2</sub> CrO <sub>4</sub>
2.	20	0				

#### Calculation of the normality of water sample (Chloride ion)

Volume of unknown water sample (V<sub>1</sub>) = 20 ml

Strength of unknown water sample (N<sub>1</sub>) = ?

Volume of standardized AgNO<sub>3</sub> (V<sub>2</sub>) = \_\_\_\_\_ ml

Strength of standardized AgNO<sub>3</sub> (N<sub>2</sub>) = \_\_\_\_\_ N

According to volumetric formula, V<sub>1</sub>N<sub>1</sub> = V<sub>2</sub>N<sub>2</sub>

$$N_1 = V_2 \times N_2 / 20$$

Strength of water sample (N<sub>1</sub>) = \_\_\_\_\_ N

#### Calculation of amount of chloride

Amount of chloride ion present in 100 ml given water sample = Eq. wt. of chloride ion x Normality of chloride ion x 100/1000 g  
 = 35.46 x Normality of chloride ion x 100/1000 g  
 = .....g

Department of Chemistry, SRMVEC

**Titration 2: Estimation of Chloride ion**

The given water sample is made up into 100 ml in a standard flask using distilled water. 20ml of this solution is pipetted out into a clean conical flask and 1ml of 2%  $K_2CrO_4$  indicator solution is added. It is then titrated against standard  $AgNO_3$  solution taken in the burette. The end point is the change of colour from yellow to reddish brown. The titration is repeated for concordant values.

**RESULT**

The amount of chloride ion present in 100 ml of the given water sample =..... g.

**Titration 1: Determination of Total Hardness**

**Sample hard water Vs Standard EDTA**

Sl.No	Volume of hard water sample (ml)	Burette reading (ml)		Volume of EDTA (ml)	Concordant value (ml)	Indicator
		Initial	Final			
1.	20	0				Eriochrome Black-T
2.	20	0				

**Calculation of the total hardness of water sample**

Volume of unknown hard water sample  $(V_1) = 20 \text{ ml}$

Strength of unknown hard water sample  $(N_1) = ?$

Volume of standard EDTA  $(V_2) = \text{_____ml}$

Strength of standard EDTA  $(N_2) = 0.01 \text{ N}$

According to volumetric formula,  
 $V_1N_1 = V_2N_2$   
 $N_1 = V_2 \times N_2 / 20$   
 $=$

Strength of hard water sample  $(N_1) = \text{_____N}$

Total hardness of the hard water sample  $= N_1 \times 50 \times 1000$

$= \text{----- ppm}$

## 2. DETERMINATION OF TOTAL, TEMPORARY & PERMANENT HARDNESS OF WATER BY EDTA METHOD

Expt. No.	
-----------	--

Date			
------	--	--	--

### AIM

To estimate the amount of total hardness, permanent hardness and temporary hardness of a given hard water sample by EDTA method. You are provided with a standard solution of calcium of strength 0.01 N.

### CHEMICALS REQUIRED

EDTA, Eriochrome black-T, water sample, ammoniacal buffer solution.

### PRINCIPLE

The amount of calcium and magnesium ions present in a given sample of water is estimated by titration with standard EDTA. The hardness is expressed in terms of calcium carbonate equivalence.

### Hardness

It is one of the important characteristics of water. Hardness of water is due to the presence of calcium and magnesium salts. Lather is not produced until the cations,  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  are completely precipitated out in the form of insoluble salts. For example,



### Temporary hardness

It is due to bicarbonates of calcium and magnesium. It can be removed by boiling. The salts are precipitated as carbonates and hydroxides respectively, which can be removed by filtration leaving only permanent hardness producing salts.

### Permanent hardness

It is due to the sulphates and chlorides of calcium and magnesium.

Department of Chemistry, SRMVEC

## Titration 2: Determination of permanent hardness

### Boiled hard water vs Standard EDTA

Sl.No	Volume of hard water (ml)	Burette reading (ml)		Volume of EDTA (ml)	Concordant value (ml)	Indicator
		Initial	Final			
1.	20	0				Eriochrome Black-T
2.	20	0				

### Calculation of the permanent hardness of water sample

Volume of unknown hard water sample  $(V_1) = 20 \text{ ml}$

Strength of unknown hard water sample  $(N_1) = ?$

Volume of standard EDTA  $(V_2) = \text{_____ml}$

Strength of standard EDTA  $(N_2) = 0.01 \text{ N}$

According to volumetric formula,  $V_1N_1 = V_2N_2$

$$N_1 = V_2 \times N_2 / 20$$

=

Strength of hard water sample  $(N_1) = \text{_____N}$

Permanent hardness of the hard water sample  $= N_1 \times 50 \times 1000$

= ----- ppm

Temporary hardness = Total hardness - Permanent hardness

Department of Chemistry, ppm/VEC


### Hardness is expressed in terms of CaCO<sub>3</sub> equivalence

The reason for choosing CaCO<sub>3</sub> as the standard for calculating hardness of water is due to its molecular weight being 100, which makes the calculation easier. It is the most insoluble salt, thus can be easily precipitated in water treatment processes.

### Determination of hardness


The hardness of water can be determined by complexometric titration. EDTA is used as a complexing agent. The Ca<sup>2+</sup> and Mg<sup>2+</sup> present in water are titrated with EDTA using Erichrome Black T as indicator.

The Chemical structure of EDTA (Ethylene diamine tetra acetic acid) is represented as below:




The indicator Eriochrome Black-T (EBT), forms a wine red colored unstable complex with calcium and magnesium ions in hard water at a pH of 8 – 10.

### STRUCTURE OF EBT



As this solution is titrated against EDTA, the Ca<sup>2+</sup> and Mg<sup>2+</sup> ions from the indicator complex form stable metal ion - EDTA complex leaving the indicator free. The free form of indicator is steel blue in colour.

The metal – EDTA complex is stable at pH 8-10. This pH range can be maintained by using ammoniacal buffer (NH<sub>4</sub>Cl + NH<sub>4</sub>OH).



The amount of EDTA used corresponds to the hardness of water.

## **PROCEDURE:**

### **Titration 1: Estimation of total hardness of sample water**

20 ml of the given hard water sample is pipetted out into a clean conical flask. 5ml of ammoniacal buffer solution and 2 drops of EBT indicator are added. The solution turns wine red in colour. The solution is then titrated against EDTA solution taken in the burette. The change of colour from wine red to steel blue is the end point. The titration is repeated for concordant values. Let the titre value be V ml. This gives the total hardness of the sample water.

### **Titration 2: Estimation of permanent hardness**

100 ml of the given sample water is taken in a clean 250 ml beaker and boiled for 10-15 minutes until it is reduced to 1/5<sup>th</sup> of the original volume. It is then cooled and filtered. The filtrate is collected in a 100 ml standard flask and made up to the mark using distilled water. 20 ml of this made up solution is pipetted out into a clean conical flask and the titration is repeated in the same way as in titration 1. Let the titre value be V ml. This gives the permanent hardness of the sample water.

### **Estimation of temporary hardness**

Temporary hardness of the water sample is calculated by subtracting permanent hardness from total hardness.

$$\text{Temporary hardness} = \text{Total hardness} - \text{Permanent hardness.}$$

## **RESULT**

1. The total hardness of the given water sample = .....ppm
2. The permanent hardness of the given water sample = .....ppm
3. The temporary hardness of the given water sample = .....ppm



**Table. 1 Titre values and different alkalinities**

S.No.	Value of [P] & [M]	Alkalinity due to			Nature of alkalinity
		OH <sup>-</sup>	CO <sub>3</sub> <sup>2-</sup>	HCO <sub>3</sub> <sup>-</sup>	
1	P = 0	0	0	M	Only HCO <sub>3</sub> <sup>-</sup> ions present
2	P = ½M	0	2P	0	Only CO <sub>3</sub> <sup>2-</sup> ions present
3	P = M	P = M	0	0	Only OH <sup>-</sup> ions present
4	P < ½ M	0	2P	M-2P	HCO <sub>3</sub> <sup>-</sup> , CO <sub>3</sub> <sup>2-</sup> ions present
5	P > ½ M	2P-M	2M – 2P	0	Only OH <sup>-</sup> , CO <sub>3</sub> <sup>2-</sup> ions present

**Table. 2 Water sample 1 vs Standard acid**

Sl.No	Volume of water sample (ml)	Burette reading (ml)			Volume of the acid (ml)	
		Initial	[P]	[M]	[P]	[M]
1.	20	0				
2.	20	0				

**Calculation**

If the data satisfies the condition  $P > 1/2M$

- (i) Volume of HCl required for OH<sup>-</sup> Alkalinity =  $2P - M = \dots\dots\dots$ ml
- (ii) Volume of HCl required for CO<sub>3</sub><sup>2-</sup> Alkalinity =  $2M - 2P = \dots\dots\dots$ ml
- (iii) HCO<sub>3</sub><sup>-</sup> absent

**Calculation of OH<sup>-</sup> alkalinity**

Volume of standardized HCl (V1) =  $\dots\dots$ ml  
 Strength of standardized HCl (N1) = 0.1 N  
 Volume of unknown water sample (V2) = 20 ml  
 Strength of unknown water sample (N2) = ?  
 According to volumetric formula,  $V_1N_1 = V_2N_2$   
 $N_2 = \dots\dots\dots V_1\text{ml} \times 0.1\text{N} / 20$   
 $= \dots\dots\dots \text{N}$

### 3. DETERMINATION OF TYPES AND AMOUNT OF ALKALINITY IN WATER SAMPLES

Expt. No.	
-----------	--

Date			
------	--	--	--

#### AIM

To determine the type and amount of alkalinity present in the given water sample. You are provided with standard solution of hydrochloric acid of strength of 0.1N.

#### CHEMICALS REQUIRED

0.1N HCL, phenolphthalein and methyl orange.


#### PRINCIPLE

The alkalinity of water is normally due to the presence of bicarbonates, carbonates and hydroxides of sodium, potassium, calcium and magnesium. The knowledge of alkalinity is important in calculating the amounts of reagents required in water softening and in internal conditioning of boiler feed waters.

The alkalinity in water can be present due to the presence of the above mentioned ions present alone or in combination as



- Hydroxide only ( $\text{OH}^-$ )
- Carbonates only ( $\text{CO}_3^{2-}$ )
- Bicarbonates only ( $\text{HCO}_3^-$ )
- Hydroxide and carbonates ( $\text{OH}^- + \text{CO}_3^{2-}$ )
- Carbonate and bicarbonates ( $\text{CO}_3^{2-} + \text{HCO}_3^-$ )

The possibility of all the three ions existing together is ruled out because the hydroxides and bicarbonate ions combine together to form  $\text{CO}_3^{2-}$  ions



Department of Chemistry, SRMVEC

The alkalinity of the water sample is determined by titrating against the standard acid using phenolphthalein and methyl orange indicators. The following reactions take place.



Amount of OH<sup>-</sup> present in 1 litre of water sample  
in terms of CaCO<sub>3</sub> equivalent = strength of water sample x Eq. wt. of CaCO<sub>3</sub>  
i.e OH<sup>-</sup> alkalinity in terms of CaCO<sub>3</sub> equivalent = .....N x 50 x 10<sup>3</sup> ppm  
Alkalinity due to OH<sup>-</sup> ion = .....ppm

**Calculation of CO<sub>3</sub><sup>2-</sup> alkalinity**

Volume of standardized HCl (V1) = .....ml  
Strength of standardized HCl (N1) = 0.1 N  
Volume of unknown water sample (V2) = 20 ml  
Strength of unknown water sample (N2) = ?  
According to volumetric formula, V<sub>1</sub>N<sub>1</sub> = V<sub>2</sub>N<sub>2</sub>  
N<sub>2</sub> = ..... V<sub>1</sub>ml x 0.1N / 20  
= .....N

Amount of CO<sub>3</sub><sup>2-</sup> present in 1 litre of water sample  
in terms of CaCO<sub>3</sub> equivalent = strength of water sample x Eq. wt. of CaCO<sub>3</sub>  
i.e CO<sub>3</sub><sup>2-</sup> alkalinity in terms of CaCO<sub>3</sub> equivalent = .....N x 50 x 10<sup>3</sup> ppm  
Alkalinity due to CO<sub>3</sub><sup>2-</sup> ion = .....ppm

**Table. 3 Water sample 2 vs Standard acid**

Sl.No	Volume of water sample (ml)	Burette reading (ml)			Volume of the acid (ml)	
		Initial	[P]	[M]	[P]	[M]
1.	20	0				
2.	20	0				

**Calculation**

Department of Chemistry, SRMVEC

If the data satisfies the condition P < 1/2M

- (i) Volume of HCl required for CO<sub>3</sub><sup>2-</sup> Alkalinity = 2P = .....ml
- (ii) Volume of HCl required for HCO<sub>3</sub><sup>-</sup> Alkalinity = M - 2P = .....ml
- (iii) OH<sup>-</sup> ion absent

The volume of acid used upto phenolphthalein end point corresponds to the reaction (1) and (2), i.e, complete neutralization of  $\text{OH}^-$  ions and neutralization of  $\text{CO}_3^{2-}$  ion up to  $\text{HCO}_3^-$  stage.

The total volume of the acid used upto methyl orange end point corresponds to all the reactions (1), (2), and (3), i.e a complete neutralization of  $\text{OH}^-$ ,  $\text{CO}_3^{2-}$  and  $\text{HCO}_3^-$  ions.

Thus from the respective volumes of acid used the strength of various ions can be determined.

## **PROCEDURE**

### **For water sample 1**

#### **Titration I: Phenolphthalein indicator**

The burette is washed, rinsed and filled with the given std. hydrochloric acid. 20 ml of the water sample (1) is pipetted into a clean conical flask. 2-3 drops of phenolphthalein indicator is added. Pink color is observed. This is titrated against the standard acid until the pink color disappears. This titre value is noted which corresponds to the phenolphthalein end point (P).

#### **Titration II: Methyl orange indicator**

Few drops of methyl orange indicator is added to the same solution after the phenolphthalein end point. The solution turns yellow in color. The titration is continued until the yellow color changes to reddish orange color. The titre value corresponds to methyl orange end point (M). The titration is repeated for concordant values. Using the titre values the type and amount of alkalinity is calculated.

### **For water sample 2**

The given water sample (2) is make up in a 100 ml standard flask with distilled water. 20ml of the made up solution is pipetted out in to a clean conical flask. The titration is repeated as above by adding phenolphthalein indicator followed by methyl orange indicator. From the titre value the alkalinity is found out.

### Calculation of $\text{CO}_3^{2-}$ alkalinity

Volume of standardized HCl	(V1)	=	.....ml
Strength of standardized HCl	(N1)	=	0.1 N
Volume of unknown water sample	(V2)	=	20 ml
Strength of unknown water sample	(N2)	=	?
According to volumetric formula,	$V_1N_1$	=	$V_2N_2$
	$N_2$	=	..... $V_1\text{ml} \times 0.1\text{N} / 20$
		=	.....N

Amount of  $\text{CO}_3^{2-}$  present in 1 litre of water sample

in terms of  $\text{CaCO}_3$  equivalent

$$= \text{strength of water sample} \times \text{Eq. wt. of } \text{CaCO}_3$$

i.e  $\text{CO}_3^{2-}$  alkalinity in terms of  $\text{CaCO}_3$  equivalent

$$= \text{.....N} \times 50 \times 10^3 \text{ ppm}$$

Alkalinity due to  $\text{CO}_3^{2-}$  ion

$$= \text{.....ppm}$$

### Calculation of $\text{HCO}_3^-$ alkalinity

Volume of standardized HCl	(V1)	=	.....ml
Strength of standardized HCl	(N1)	=	0.1 N
Volume of unknown water sample	(V2)	=	20 ml
Strength of unknown water sample	(N2)	=	?
According to volumetric formula,	$V_1N_1$	=	$V_2N_2$
	$N_2$	=	..... $V_1\text{ml} \times 0.1\text{N} / 20$
		=	.....N

Amount of  $\text{HCO}_3^-$  present in 1 litre of water sample

in terms of  $\text{CaCO}_3$  equivalent

$$= \text{strength of water sample} \times \text{Eq. wt. of } \text{CaCO}_3$$

i.e  $\text{HCO}_3^-$  alkalinity in terms of  $\text{CaCO}_3$  equivalent

$$= \text{.....N} \times 50 \times 10^3 \text{ ppm}$$

Alkalinity due to  $\text{HCO}_3^-$  ion

$$= \text{.....ppm}$$

Department of Chemistry, SRMVEC

## RESULT

The given water sample (1) contains

1. Hydroxide alkalinity = .....ppm
2. Carbonate alkalinity = .....ppm
3. Bicarbonate alkalinity = .....ppm


The given water sample (2) contains

1. Carbonate alkalinity = .....ppm
2. Bicarbonate alkalinity = .....ppm
3. Hydroxide alkalinity = .....ppm

### Titration of mixture of BaCl<sub>2</sub> Vs Na<sub>2</sub>SO<sub>4</sub>

S.No	Volume of Na <sub>2</sub> SO <sub>4</sub> added (ml)	Conductance (mho)
1.		
2.		
3.		
.		
.		
.		
20.		

#### Model Graph



#### Calculation of strength of BaCl<sub>2</sub>

Volume of unknown BaCl<sub>2</sub> (V<sub>1</sub>) = 20 ml

Strength of unknown BaCl<sub>2</sub> (N<sub>1</sub>) = .....N

Volume of standard Na<sub>2</sub>SO<sub>4</sub>(V<sub>2</sub>) = ..... ml (titre value)

Strength of standard Na<sub>2</sub>SO<sub>4</sub>(N<sub>2</sub>) = 0.1 N

Strength of unknown BaCl<sub>2</sub> (N<sub>1</sub>) = V<sub>2</sub>N<sub>2</sub> / V<sub>1</sub>

Department of Chemistry, SRMVEC

Strength of unknown BaCl<sub>2</sub> (N<sub>1</sub>) = .....N

The amount of barium chloride present in 1 litre of the given solution

= strength of BaCl<sub>2</sub> x Equivalent weight of BaCl<sub>2</sub>

= ..... N x 122.14

= ..... g/L

#### 4. CONDUCTOMETRIC TITRATION OF BARIUM CHLORIDE AGAINST SODIUM SULFATE (PRECIPITATION TITRATION)

Expt. No.	
-----------	--

Date			
------	--	--	--

##### AIM


To estimate the amount of BaCl<sub>2</sub> present in the given solution. You are provided with standard solution of Na<sub>2</sub>SO<sub>4</sub> of strength 0.1 N.

##### MATERIALS REQUIRED

Conductivity meter, conductivity cell, micro burette, pipette, standard 0.5 N Na<sub>2</sub>SO<sub>4</sub>, BaCl<sub>2</sub> solution.

##### PRINCIPLE

Conductivity of the solution depends upon number of ions in solution. During the addition of Na<sub>2</sub>SO<sub>4</sub> to BaCl<sub>2</sub> taken in the beaker. Ba<sup>2+</sup> ions are removed as BaSO<sub>4</sub> precipitate (which produce negligible amount of ions).



Hence the conductance decreases sharply. When all the Ba<sup>2+</sup> ions are removed as BaSO<sub>4</sub>, further addition of Na<sub>2</sub>SO<sub>4</sub> produces 2Na<sup>+</sup> and SO<sub>4</sub><sup>2-</sup> ions. Hence conductivity starts increasing sharply. The intersection point corresponds to the end point of the titration.

##### PROCEDURE

Sodium sulphate is taken in the burette. 20ml of the given BaCl<sub>2</sub> solution is pipetted out into a clean beaker. Conductivity cell is placed in it and diluted with conductivity water.

1 ml of Na<sub>2</sub>SO<sub>4</sub> solution is added to the solution, stirred and then conductance is measured. This is continued up to the end point. The conductance decreases up to the end point. After the end point, conductance increases due to the addition of excess of Na<sub>2</sub>SO<sub>4</sub>.

A graph is plotted between the volume of Na<sub>2</sub>SO<sub>4</sub> and conductance. From the graph end point is noted and the amount of BaCl<sub>2</sub> present in 1 litre of the solution is calculated.

##### RESULT

The amount of BaCl<sub>2</sub> present in 1 litre of the given solution = ..... g/L



### Titration 1: Determination of acidity of water sample 1

#### Sample water Vs Standard Na<sub>2</sub>CO<sub>3</sub>

Sl.No	Volume of water sample 1 (ml)	Burette reading (ml)		Volume of Standard Na <sub>2</sub> CO <sub>3</sub> (ml)	Concordant value (ml)	Indicator
		Initial	Final			
1.	20	0				Methyl orange
2.	20	0				

#### Calculation of acidity of water sample

Volume of water sample (V<sub>1</sub>) = 20 ml

Strength of water sample (N<sub>1</sub>) = ?

Volume of standard Na<sub>2</sub>CO<sub>3</sub> (V<sub>2</sub>) = \_\_\_\_\_ ml

Strength of standard Na<sub>2</sub>CO<sub>3</sub> (N<sub>2</sub>) = 0.01 N

According to volumetric formula, V<sub>1</sub>N<sub>1</sub> = V<sub>2</sub>N<sub>2</sub>

$$N_1 = V_2 \times N_2 / 20$$

Strength of water sample 1 (N<sub>1</sub>) = \_\_\_\_\_ N

Amount of acidity present in 1 litre of water sample in

terms of CaCO<sub>3</sub> equivalent

$$= \text{strength of water sample} \times \text{Eq. wt. of CaCO}_3 \times 1000$$

Department of Chemistry, SRMVEC

$$= \dots\dots\dots N \times 50 \times 10^3 \text{ ppm}$$

The acidity of given water sample 1

$$= \dots\dots\dots \text{ppm}$$

## 5. PREPARATION OF Na<sub>2</sub>CO<sub>3</sub> AS A PRIMARY STANDARD AND ESTIMATION OF ACIDITY OF A WATER SAMPLE USING THE PRIMARY STANDARD

Expt. No.	
-----------	--

Date			
------	--	--	--

### AIM

To determine the acidity of water sample.

### MATERIALS REQUIRED

Apparatus:

100 mL beaker, 100 mL standard flask, 250 mL conical flask, 20 mL pipette, burette, retort stand

Chemicals:



Na<sub>2</sub>CO<sub>3</sub>, concentrated HCl, Methyl orange indicator

### PRINCIPLE

Acidity of a water is its quantitative capacity to react with a strong base to a designated pH. The measured value may vary significantly with the end-point pH used in the determination. Acidity is a measure of aggregate property of water and can be interpreted in terms of specific substances only when the chemical composition of the sample is known. Strong minerals acids, weak acids such as carbonic and acetic, and hydrolyzing salts such as iron or aluminium sulfate may contribute to the measured acidity according to the method of determination.

Pure sodium carbonate is a non-hygroscopic primary standard reagent. A standard solution prepared by dissolving a known amount of the solid Na<sub>2</sub>CO<sub>3</sub> in a fixed volume can be used to standardize other acid solutions. The aim of the experiment is to prepare a standard solution of exact concentration which can later be used to standardize other solutions.

HCl is not a primary standard. Thus, after a dilute HCl solution is prepared, it has to be standardized with a primary standard solution in order to determine its concentration accurately. In this experiment a standard Na<sub>2</sub>CO<sub>3</sub> solution and a dilute HCl solution will be prepared. The HCl solution will later be standardized against the Na<sub>2</sub>CO<sub>3</sub> solution. The reaction between HCl and Na<sub>2</sub>CO<sub>3</sub> is as follows.



These are the hydrochloric acid and sodium carbonate equations. In acid-base titrations, the amount of the acid becomes chemically equivalent to the amount of base present. The solution becomes neutral in the event of a strong acid and a strong base titration of the solution.

**Titration 2: Determination of acidity of water sample 2**

**Sample water Vs Standard Na<sub>2</sub>CO<sub>3</sub>**

Sl.No	Volume of water sample 2 (ml)	Burette reading (ml)		Volume of Standard Na <sub>2</sub> CO <sub>3</sub> (ml)	Concordant value (ml)	Indicator
		Initial	Final			
1.	20	0				Methyl orange
2.	20	0				

**Calculation of acidity of water sample**

Volume of unknown water sample (V<sub>1</sub>) = 20 ml  
 Strength of unknown water sample (N<sub>1</sub>) = ?  
 Volume of standard Na<sub>2</sub>CO<sub>3</sub> (V<sub>2</sub>) = .....ml  
 Strength of standard Na<sub>2</sub>CO<sub>3</sub> (N<sub>2</sub>) = 0.01 N  
 According to volumetric formula, V<sub>1</sub>N<sub>1</sub> = V<sub>2</sub>N<sub>2</sub>  


$$N_1 = V_2 \times N_2 / 20$$
 Strength of water sample 2 (N<sub>1</sub>) = .....N

Amount of acidity present in 1 litre of water sample in terms of CaCO<sub>3</sub> equivalent = strength of water sample x Eq. wt. of CaCO<sub>3</sub> x 1000


Department of Chemistry, SRMISTC = .....N x 50 x 10<sup>3</sup> ppm

The acidity of given water sample 2 = .....ppm

The chemical structure of methyl orange is represented as below:



The indicator **methyl orange** (MeO), forms a red colored complex with acidic water sample a pH of 3-5.



As this solution is titrated against  $\text{Na}_2\text{CO}_3$ , the indicator forms a stable complex leaving the indicator free. The free form of indicator is red in colour.

## PROCEDURE

### Preparation of standard $\text{Na}_2\text{CO}_3$ solution

1. Weigh accurately about 1.33 g of  $\text{Na}_2\text{CO}_3$ . Record the exact mass of the  $\text{Na}_2\text{CO}_3$ .
2. Dissolve the  $\text{Na}_2\text{CO}_3$  in about 50 mL water in a 100 mL beaker.
3. Dissolve the  $\text{Na}_2\text{CO}_3$  in about 50 mL water in a 100 mL beaker.
4. Transfer the solution into a 250 mL volumetric flask. Rinse the beaker with distilled water.
5. Add distilled water to the mark. Stopper the flask. Shake it by turning it upside down a few times to make the solution homogeneous.

Department of Chemistry, SRMVEC

### Titration: Determination of acidity of water sample 1

The burette is washed well with distilled water and rinsed with small amount of  $\text{Na}_2\text{CO}_3$  solution. The pipette is washed with distilled water and rinsed with small amount of water sample. 20ml of water sample is pipetted out into a clean conical flask. 1 ml or 3 drops of methyl orange indicator is added and titrated against standard  $\text{Na}_2\text{CO}_3$  solution taken in the burette. The end point is the change of colour from red orange to yellow. The titration is repeated for concordant values. From the titre value, the acidity is found out.

**Titration: Determination of acidity of water sample 2**

The given water sample is made up into 100 ml in a standard flask using distilled water. 20ml of the made up water sample is pipetted out in to a clean conical flask. This is titrated against the standard  $\text{Na}_2\text{CO}_3$  until the indicator colour changes from red orange to yellow. The titration is repeated for concordant values. From the titre value, the acidity is found out.

**RESULT**

1. The acidity of given water sample 1 = .....ppm
2. The acidity of given water sample 2 = .....ppm

**CALCULATION:**

Subtract the initial weight (in grams) of the empty container from the weight of the container with the dried residue to obtain the increase in weight.

A- Weight of clean dried container (g)

B- Weight of container and residue (g)

C- Volume of Sample (ml)

Sl.No	Volume of water sample (ml)	Weight of clean dried container (g)	Weight of container and residue (g)	Weight of residue solids (mg/g)
1.				
2.				

Dissolved solids, TDS (mg/l)

= mg of solids in the beaker / (volume of sample)x 1000

= ((B - A) / C) \* (1000 mg/g) / (1000) ml/L

= \_\_\_\_\_ mg/L

## 6. ESTIMATION OF TDS OF A WATER SAMPLE BY GRAVIMETRY

Expt. No.	
-----------	--

Date			
------	--	--	--

### AIM

To determine the Total dissolved solids (TDS) of given water samples.

### MATERIALS REQUIRED

Water sample, electronic balance machine, beaker, measuring cylinder, Whatman filter paper/Gooch crucible or sintered crucible, funnel, dropper, evaporating dish/ceramic dish, desiccator

### PRINCIPLE

#### *Introduction:*

Water is a good solvent and picks up impurities easily. Pure water - tasteless, colorless, and odorless is often called the universal solvent. Dissolved solids" refer to any minerals, salts, metals, cations or anions dissolved in water. Total dissolved solids (TDS) comprise inorganic salts (principally calcium, magnesium, potassium, sodium, bicarbonates, chlorides, and sulfates) and some small amounts of organic matter that are dissolved in water.

In general, the total dissolved solids concentration is the sum of the cations (positively charged) and anions (negatively charged) ions in the water. Therefore, the total dissolved solids test provides a qualitative measure of the amount of dissolved ions but does not tell us the nature or ion relationships

An elevated total dissolved solids (TDS) concentration is not a health hazard. The TDS concentration is a secondary drinking water standard and, therefore, is regulated because it is more of an aesthetic rather than a health hazard. An elevated TDS indicates the following:

- 1) The concentration of the dissolved ions may cause the water to be corrosive, salty or brackish taste, result in scale formation, and interfere and decrease efficiency of hot water heaters; and
- 2) Many contain elevated levels of ions that are above the Primary or Secondary Drinking Water Standards, such as an elevated level of nitrate, arsenic, aluminum, copper, lead, etc.

Department of Chemistry, SRMVEC

### PROCEDURE

1. Filter your water sample through a Whatman Filter paper.
2. Collect the filtrate (liquid) and rinse water in a beaker/flask.
3. Take the weight of empty container (beaker/ ceramic dish/ evaporating dish). Make sure the container should be dried.
4. Add the filtrate to the container and allow the sample to stay in the oven at 105°C for 24 hours. If

possible, increase the temperature of the drying oven to 180°C and allow the sample to dry for up to 8 hours.

5. Remove the container - Remember it is very hot. After removing from the drying oven, the sample container should be placed in a desiccator to cool in a dry air environment for at least 3 to 4 hours.
6. After the container cools, reweigh the container at least three times. Then do the following:

## **RESULT**

The dissolved solids, TDS of given water sample = ..... (mg/L) or ppm



## 7. PREPARATION OF NANOPARTICLES (AgNPs)

Expt. No.	
-----------	--

Date			
------	--	--	--

### AIM


To synthesize the silver nanoparticles (AgNPs) by chemical reduction method..

### MATERIALS REQUIRED

Silver nitrate (AgNO<sub>3</sub>), serine, polyvinyl pyrrolidone (PVP), and NaOH.

### PRINCIPLE

The Ag NPs were synthesized by the reduction of Ag<sup>+</sup> ions in the presence of microwave radiation using serine (amino acid) as a reducing agent and PVP as a capping agent. There is an oxidation of -OH group of serine into aldehyde group by reducing Ag<sup>+</sup> and finally leading to the formation of Ag<sup>0</sup> nuclei, which can be stabilized by PVP.



### PROCEDURE

Synthesis of silver nanoparticles

A microwave oven was used in the typical synthesis of Ag NPs. 100 mL of 5 mM silver nitrate solution was prepared in 250 mL conical flask using the de-ionized water. 0.0525 g of serine and 0.75 g of PVP were added into the above silver nitrate solution, and then 4 drops of 1 M NaOH were added by altering the pH for the favorable reduction. The above mixed solution was kept in the microwave oven (input power 1200 W, 50 Hz) for 90 s. The solution turned into yellowish brown colour, which confirmed the formation of Ag NPs.

Department of Chemistry, SRMVEC

### RESULT

The AgNPs were successfully synthesized by chemical reduction method with help of microwave oven.